LOUISIANA DEPARTMENT OF TRANSPORTATION & DEVELOPMENT VERIFICATION PROCEDURE

DRYING OVEN

Equipment

- 1. A calibrated thermometer in 1.0 °C increments having a range which includes the temperature range to be checked.
- 2. A brass thermometer well.
- 3. A clothes pin on the thermometer.

Procedures

- 1. Place the clothes pin on the thermometer.
- 2. Place the thermometer inside the brass well with the clothes pin attached. Position the thermometer, brass well and clothes pin on the shelf where the samples are normally dried.
- 3. Take the first reading at least 1 hour after closing the oven (oven should remain undisturbed).
- 4. Take three consecutive readings no less than 1/2 hour apart to determine if the temperature range is within the specified tolerance.
- 5. Adjust the temperature of the oven if an observed temperature reading is outside the tolerance specified (allow at least 1/2 hour for the temperature to stabilize between each adjustment). Return to step C.

STATE OF LOUISIANA DEPARTMENT OF TRANSPORTATION AND DEVELOPMENT

VERIFICATION OF OVEN TEMPERATURE SETTINGS verification procedure used: DOTD A1

| Verification frequency: 4 months | Previous verification date: | |
|-------------------------------------|------------------------------|--|
| Date of verification: | Next verification due: | |
| Identification no.: | Mfg./distributor: | |
| Verified by: | Verification equipment used: | |
| Reference to certified thermometer: | YesNo | |

| Oven temperature setting | Oven temperature readout | Scale |
|--------------------------|--------------------------|-------|
| | | |

| Time of verification | Temperature Readings | Scale | Remarks |
|----------------------|-------------------------|-------|---------|
| | | | |
| | | | |

Recommended action: Repair_____ Replace _____ None _____

Comments: